Application No.: 10/632,039 Attorney Docket No.: JJM-346CIP1 Response to Office Action dated January 26, 2009

## CLAIMS

The Listing of Claims replaces all prior versions of claims in the subject application.

 (Currently amended) A method for cleaning and sterilizing a medical device comprising:

placing the device into a container:

cleaning the device in the container with a cleaning solution;

rinsing the device in the container with a rinse solution;

vaporizing a liquid substance in the container to create a sterilant vapor, wherein the liquid substance comprises a retained portion of the rinse solution; and

contacting the device with the vapor to effect sterilization of the device.

- 2. (Previously presented) The method according to claim 1 further comprising storing the device in the container in sterile form.
- 3. (Canceled)
- (Previously presented) The method according to claim 1 wherein the rinse solution comprises a chemical sterilant.
- (Previously presented) The method according to claim 4 wherein the chemical sterilant comprises hydrogen peroxide.

Application No.: 10/632,039 Attorney Docket No.:JJM-346CIP1 Response to Office Action dated January 26, 2009

(Previously presented) The method according to claim 1 wherein the liquid substance comprises a chemical sterilant.

- (Previously presented) The method according to claim 6 wherein the chemical sterilant comprises hydrogen peroxide.
- (Previously presented) The method according to claim 1 and further comprising introducing the liquid substance as a mist.
- (Previously presented) The method of claim 1 further comprising retaining a
  predetermined amount of the liquid substance in the container prior to the
  vaporizing.
- 10. (Currently amended) A method for cleaning and sterilizing a medical device comprising:

placing the device into a container;

device

cleaning the device in the container with a cleaning solution;

rinsing the device in the container with a rinse solution;

retaining a predetermined amount of a liquid substance in the container, wherein the liquid substance comprises a retained portion of the rinse solution;

vaporizing the liquid substance in the container to create a vapor; and contacting the device with the vapor to effect sterilization of the

Application No.: 10/632,039
Attorney Docket No.:JJM-346CIP1
Response to Office Action dated January 26, 2009

 (Previously presented) The method according to claim 10 wherein the rinse solution comprises a chemical sterilant.

- 12. (Previously presented) The method according to claim 11 wherein the chemical sterilant comprises hydrogen peroxide.
- (Previously presented) The method according to claim 10 wherein the liquid substance comprises a chemical sterilant.
- 14. (Previously presented) The method according to claim 13 wherein the chemical sterilant comprises hydrogen peroxide.
- 15. (Previously presented) A method for cleaning and sterilizing a medical device comprising:

placing the device into a container;

cleaning the device in the container with a cleaning solution;

rinsing the device in the container with a rinse solution;

retaining a portion of a liquid substance in the container, wherein the liquid substance comprises at least a portion of the rinse solution;

vaporizing the liquid substance in the container to create a sterilant vapor; and

contacting the device with the vapor to effect sterilization of the device.

16. (Canceled)

Application No.: 10/632,039 Attorney Docket No.:JJM-346CIP1 Response to Office Action dated January 26, 2009

- 17. (Previously presented) The method according to claim 15 wherein the rinse solution comprises a chemical sterilant.
- 18. (Previously presented) The method according to claim 17 wherein the chemical sterilant of the rinse solution comprises hydrogen peroxide.
- 19. (Previously presented) The method according to claim 10, wherein at least a portion of the liquid substance comprises water.
- 20. (Previously presented) The method according to claim 15, wherein at least a portion of the liquid substance comprises water.